

Arita Meibomian

GLAND COMPRESSOR



This new instrument has been designed by Dr. Reiko Arita of Saitama City, Japan to gently express meibum from dysfunctional meibomian glands (MGD). The special angles of the forceps allow the surgeon to easily insert and position the tips over the eyelid parallel to the lid margin. The broad, smoothly polished jaws are then used to atraumatically compress the glands – relieving any occlusions with minimal pain to the patient. This instrument is user-friendly while addressing an upper or lower eyelid and working from either a temporal or nasal approach.



Watch it!

K5-5900

 **katena**
DESIGNED FOR SIGHT™